


# In-Fusion<sup>®</sup> HD Cloning Plus



Fast, accurate cloning for simple to complex projects



Clontech TakaRa cellortis

# In-Fusion<sup>®</sup> Cloning:



# Ideal for...

- Single-insert cloning
- Multiple-insert cloning
- Large-insert cloning
- Adding small tags
- Direct cloning into large vectors (up to 46-kb cosmids)
- Site-directed mutagenesis (single or multiple deletions, insertions, or substitutions)
- HTP Cloning


# The right clone every time.

Good science relies on a solid foundation of inspiration matched with technical expertise. Every decision has an impact, starting with experimental design. This is why we are dedicated to delivering best-in-class products backed by expert scientific support. We understand that streamlined methods and high-performance reagents enable you to push your research to its fullest potential.

In-Fusion Cloning kits are unique in their ability to allow easy, accurate, seamless cloning without the use of ligase. From the simple to the complex, you can recover your final construct on the second day! In just one quick reaction, you can clone ANY insert into ANY vector at ANY locus, without the hassle of subcloning, unwanted extra bases, or inconvenient restriction sites.

The reaction is specific and directional, ensuring an exceptionally high rate of cloning accuracy in all applications. Our standards for accuracy involve more than simple transformation efficiency. Why waste time screening a large number of colonies? With background levels approaching zero, you can pick just a few colonies and be confident that you have the right clone.

As the first step in so many experiments, successful cloning is key to recovering meaningful data. The cloning tools you use can be the difference between wasted time and fast, consistent results.






#### Figure 1. In-Fusion Cloning protocol.

The key to In-Fusion Cloning is 15 bp of homology between your insert(s) and linearized vector backbone. Add these 15 bp to the 5' end of your forward and reverse primers, and PCR amplify your desired insert(s). Linearize your vector, by restriction digest or inverse PCR. Combine in a single tube containing the In-Fusion enzyme mix, perform the In-Fusion cloning reaction, and transform into competent cells.

# One Cloning System for multiple applications

### Clone multiple fragments with ease!

Using In-Fusion, cloning multiple fragments simultaneously is as easy as cloning a single fragment and it can be done in a single reaction! Just combine the PCR fragments with appropriate complementary ends and a linear vector, then incubate with the In-Fusion enzyme. This saves weeks and months that would otherwise be spent screening clones and sub-cloning. The ability to easily, rapidly and precisely clone many fragments at once will help to speed up the generation of complex target constructs in your lab. And don't just take our word for it; Zhu et al. successfully constructed a three domain immunoglobulin fusion protein by means of a four-piece In-Fusion reaction (See Fig. 2 below adapted from Zhu *et al.* (2007) BioTechniques 43(3): 354-359.).




# Figure 2 : Seamless construction of a three-domain immunoglobulin fusoin protein with a four piece In-Fusion reaction.

Segments were generated by PCR using primers that contained a 15 bp overlap with the adjacent segment, and 20-30 bp of segment specific sequence. Colored regions indicate overlap regions with 15 bp of identify. Arrows indicate PCR primers. The segments were joined to a Ncol-Sall-digested expression vector in a ligase-free In-Fusion reaction.

#### Site-directed mutagenesis

In-Fusion cloning can be used for both cloning and site-directed mutagenesis including insertions, deletions and substitutions, eliminating the need for a separate kit for mutagenesis. Save money, time, and space in your laboratory by using this one universal system for genetic manipulation, including the engineering of both single and multiple mutations quickly and efficiently in a single reaction. Mutagenesis can be performed easily by combining the strength of the In-Fusion enzyme with inverse PCR using mutagenic primers (See Fig. 3 below).

Change occurs here



### HTP Cloning for protein expression and production studies

Obtaining recombinant proteins in a soluble form suitable for crystallization remains a bottleneck for high-throughput (HTP) structural biology. The Oxford Protein Production Facility (OPPF) has created the pOPIN vector suite for use with Clontech's In-Fusion cloning system. This custom vector suite in combination with the In-Fusion enzyme delivers a versatile, precise, and seamless one-step cloning process (i.e., no unwanted residues in the translated protein), that can be readily adapted to high-throughput structural biology studies. Incorporating the In-Fusion system into the OPPF HTP pipeline enabled parallel vector construction, expression screening, and purification, accelerating the protein production process (Berrow *et al.*, 2007, Nucleic Acids Res. 35(6):e45).

# In-Fusion HD Cloning Plus Systems

### **Complete, optimized systems**

Each system fully equips you for every step, bringing the price per reaction in line with the true total cost of your experiment. The components have been optimized to work together efficiently, and provide you with a high out-of-the-box success rate.

Our complete range of In-Fusion HD Cloning Plus kits allows you to choose the version that best fits your experimental needs.

The flexibility and speed of the In-Fusion mechanism means it can be easily adapted for nearly every cloning application.

## Complete systems include:

- CloneAmp<sup>™</sup> HiFi PCR Premix
- **PCR product treatment** (either Cloning Enhancer or a NucleoSpin Gel and PCR Cleanup kit)
- In-Fusion HD Enzyme Premix
- Stellar<sup>™</sup> Competent Cells

Be Green!

The In-Fusion HD EcoDry System can be stored at room temperature.

	In-Fusion HD	In-Fusion HD Cloning	In-Fusion HD EcoDry™
	Cloning Plus System	Plus CE System	Cloning Plus System
System format	liquid	liquid	lyophilized
Ability to customize reaction volume and/or plasticware	$\checkmark$	$\checkmark$	
Easily scalable for high-throughput applications		$\checkmark$	
Minimize handling with pre-aliquoted, lyophilized components			$\checkmark$
Store reagents at room temperature			$\checkmark$
Fast, enzymatic removal of background DNA and PCR residue (for PCR that produces a single, clean band)		$\checkmark$	
Suitable for PCR and gel extraction when PCR produces >1 product	$\checkmark$		$\checkmark$



### Exceptional cloning accuracy with In-Fusion HD Cloning.

Consider the following scenario:

#### **The Experiment**

One vector Three inserts One 15-minute reaction

#### The Results

26/26 colonies (100% cloning accuracy) Sequence-verified Low background (single colony on no-insert plate)

In-Fusion technology reliably demonstrates high cloning accuracy in all applications. This is especially useful with more challenging projects such as multiple-insert cloning, where high background and false positives present a bigger problem.

## Table I. Selection Guide for In-Fusion HD Cloning Systems


# PCR Polymerases

A crucial element required for successful cloning is an error-free insert. A PCR polymerase capable of amplifying your templates for both cloning and colony screening is key to seeing accurate results. Clontech offers a wide range of DNA polymerases to suit every cloning project.

- CloneAmp HiFi PCR Premix
- SapphireAmp<sup>®</sup> Fast PCR Master Mix
- PrimeSTAR® GXL DNA Polymerase


#### Table III. Selection Guide for DNA Polymerases for Cloning

CloneAmp HiFi PCR Premix	SapphireAmp Fast PCR Master Mix	PrimeSTAR GXL DNA Polymerase
Optimized for use with In-Fusion Cloning Extremely low error rate (see Fig.4 below) Quick priming & extension (5sec/kb) Templates up to 10 kb GC-rich templates	Analyze without purification Restriction digest Gel electrophoresis Colony PCR in ~1 hour Master mix includes dye (see Fig.5 below)	Robust, high-fidelity PCR in all conditions Excels with long templates GC- or AT-rich targets Templates with repeats



#### Figure 4 : Mutation frequency of CloneAmp Hifi compared to other high-fidelity PCR enzymes

Mutation frequency of CloneAmp HiFi Polymerase compared to other high-fidelity PCR enzymes. Eight arbitrarily selected GC-rich regions were amplified with CloneAmp HiFi Polymerase or other DNA polymerases using a *Thermus thermophilus* HB8 genomic DNA template, and cloned into suitable plasmids. Multiple clones were selected for each amplification product and subjected to sequence analysis. DNA fragments amplified using CloneAmp HiFi Polymerase yielded only 12 mismatched bases per 542,580 total bases—lower than an alternative high-fidelity enzyme from Company A, and 10-fold lower than Taq DNA polymerase.



#### Figure 5. Colony PCR protocol with SapphireAmp Fast PCR Master Mix.

Colony PCR can be set up with two basic steps: 1. Pick a colony ; 2. Inoculate a replica plate.

#### PRODUCTS

Cat. #	Product	Size
In-Fusion Cl	oning	
638909	In-Fusion HD Cloning Plus	10 Rxns
638910	(Liquid system, includes: In-Fusion HD Enzyme Premix, CloneAmp HiFi PCR Premix,	50 Rxns
638911	Stellar Competent Cells, and NucleoSpin Gel and PCR Clean-Up kit)	100 Rxns
638920		96 Rxns
638916	In-Fusion HD Cloning Plus CE	10 Rxns
638917	(Liquid system, includes: In-Fusion HD Enzyme Premix, CloneAmp HiFi PCR Premix,	50 Rxns
638918	Stellar Competent Cells, and Cloning Enhancer)	100 Rxns
638919		96 Rxns
638912	In-Fusion HD EcoDry Cloning Plus	8 Rxns
638913	(Lyophilized system, includes: In-Fusion HD Enzyme Premix, CloneAmp HiFi PCR Premix,	24 Rxns
638914	Stellar Competent Cells, and NucleoSpin Gel and PCR Clean-Up kit)	48 Rxns
638915		96 Rxns
PCR Polyme	rases	
639298	CloneAmp HiFi PCR Premix	200 Rxns
RR350A	SapphireAmp Fast PCR Master Mix	160 Rxns
RR350B		800 Rxns
R050A	PrimeSTAR GXL DNA Polymerase	250 Units
R050B		1,000 Units

#### Notice to Purchaser.

Your use of these products and technologies is subject to compliance with any applicable licensing requirements described on the product's web page at http://www.takarabio.com. It is your responsibility to review, understand and adhere to any restrictions imposed by such statements.

Clontech® products are to be used for research purposes only. They may not be used for any other purpose, including, but not limited to, use in drugs, *in vitro* diagnostic purposes, therapeutics, or in humans. Clontech products may not be transferred to third parties, resold, modified for resale, or used to manufacture commercial products or to provide a service to third parties without prior written approval of Takara Bio USA, Inc.



#### Takara Bio Europe

ordersEU@takarabio.com • techEU@takarabio.com • infoEU@takarabio.com • +33 139 046 880 For Research Use Only. Not for use in diagnostic or therapeutic procedures. Not for resale. Takara, the Takara logo, SapphireAmp, and PrimeSTAR are trademarks of TAKARA HOLDINGS, Kyoto, Japan. Clontech, the Clontech logo, CloneAmp, EcoDry, In-Fusion, Stellar, and that's GODD scienceI are trademarks of Takara Bio USA, Inc. All other marks are the property of their respective owners. Certain trademarks may not be registered in all jurisdictions. ©2017 Takara Bio Europe