#### TWIST BIOSCIENCE

# Human Core Exome Kit

**Twist Bioscience**, the leader in synthetic DNA with unparalleled precision at scale, is redefining targeted sequencing performance with the Twist Human Core Exome Kit. This kit includes library preparation and enrichment components, allowing the flexibility to lower sequencing costs, increase sample throughput, and achieve a higher depth of coverage across target regions with uncompromising quality.

#### EXOME DEFINITION

- 33 Mb CCDS coverage
- 99% of ClinVar variants covered
- 99.3% of targeted regions covered at 20x with 5.3 Gb

#### **KEY BENEFITS**

Unparalleled uniformity for sequencing efficiency

Lower overall cost of sequencing

# Deeper target coverage of difficult regions

High complexity libraries for high sensitivity variant calling

#### **Greater flexibility**

Seamless integration into any workflow with scalable throughput




Figure 1. Twist Human Core Exome Kit resulted in a higher number of targets covered with 20 or 30 reads compared to leading competitors with comparable coverage of 1-Plex and 8-Plex. This result was achieved with a 1.32 fold 80 base penalty (green) compared to Competitor I's (grey) 1.65. The graphs also show high on-target rate and 50% less duplicates compared to Competitor I.



## **Revolutionizing Target Enrichment**

A critical component for successful sequencing analyses is the quality of the capture reaction, which is driven by the uniformity of the reads. A more uniform capture leads to more uniform sequencing as shown by the lower fold 80 base penalty. With higher uniformity, sequencing reads can be effectively utilized to drive coverage on targets, limiting wasted reads.

High uniformity is the combined result of the probe design algorithm, dsDNA probes, and probe boosting based on GC content, amplification and enrichment bias. NGS QC on all Twist probe pools ensures that all probes are present at the correct concentrations and reproducibility of results is independent of our production.

	A-1	A-2	I-1	TWIST
	1-PLEX	1-PLEX	8-PLEX	8-PLEX
FOLD_80_BASE_PENALTY	1.93	1.71	1.69	1.34
HS_LIBRARY SIZE	123M	75M	123M	644M
PCT_EXC_DUPE	8.4%	11.3%	10.0%	1.6%





Figure 2. Twist Bioscience provides more uniformity and library complexity than competitor kits. Sequencing depth is 150x raw coverage

## **NGS Workflow**



# **Customization of Core Exome Kit Content and Components**

The ability to simply add your spike-in panel to the Core Exome creates added differentiation based on your research requirements. Rapid design optimization with a quick production turnaround reduces overall development time for your custom exome.



# Your Research. Your Needs. Your Kit.

Not all labs have the same workflow needs, so why settle for a one-size-fits-all solution? Twist Human Core Exome Kit gives you the flexibility to select components in a modular fashion. With four workflow configurations, one can utilize the complete kit or easily add components with minimal impact on an existing process.

With Twist Human Core Exome Kit, running more samples per sequencing run doesn't mean compromising on quality. When using Twist Human Core Exome Kit, single and multiplex enrichment efficiencies are very similar. Don't alter a workflow for a product, let Twist's NGS Solutions complement any research needs.




Figure 5. The versatility of a modular kit allows researchers to purchase the desired components needed for their workflow. Here, each kit contains reagents for 16 and 96 samples from singleplex to 8-plex hybridizations

# **Don't Settle for Less**

DISCOVER THE TWIST ADVANTAGE



### **HIGH QUALITY PROBE SYNTHESIS**

Unsurpassed uniformity and low error rate with Twist's silicon-based DNA synthesis platform enable high library performance metrics and singleplex data quality with multiplex capture



## IMPROVED DESIGN

Focus on relevant regions and save on sequencing cost without missing important data




## dsDNA PROBES

Target both strands of the genome for improved sensitivity



### **NGS-BASED QC ON FINALIZED PROBES SET**

Remove uncertainty and guesswork from your assays



### FLEXIBLE AND CUSTOMIZABLE WORKFLOW

Complete kit and modular library preparation and enrichment reagents for workflow integration with minimal optimization, and easy incorporation of custom content

These Products are subject to certain use restrictions as set forth in Twist's Supply Terms and Conditions www.twistbioscience.com/supply-terms-and-conditions

