

A Complete Portfolio for Your DNA/RNA Seq WORKFLOW

18 giugno 2019

CREO Perugia

Sala riunioni secondo piano, 11:00 - 12:00

Program:

Speaker: Pierluigi Tenca - Sales Manager Diatech Labline

11:00 - 12:00

RNA-seq SMARTer kits:

- Single Cell Total RNA-seq kits (from 1 to 1000 cells or 10 pg-10 ng of Total RNA, RIN 2-8)
- Stranded Total RNA-seq kits (RIN 2-10, from 200 pg -1 µg of RNA)
- Small RNA-seq kit (from 1 ng-2 ug of RNA; gel-free size selection)
- BCR RNA-seq kit

DNA-seq kits:

- SMARTer ThruPLEX DNA-seq (from 50pg to 50ng of gDNA, cDNA, FFPE DNA, ChIP DNA)
- SMARTer ThruPLEX Plasma-seq (1ng to 30ng cfDNA from plasma and urine, liquid biopsy)
- SMARTer ThruPLEX Tag-seq (UMTs for rare allele detection in cancer and other diseases);
- PicoPLEX Gold DNA-seq kit (from Single cell or 6 pg of DNA input)

Microbiome 16S&18S long Reads-Seq using Illumina Platform:

- LoopSeq 16S&18S Long Reads

2S Turbo Seq – Real Tagmentation Alternative

Targeted Exome and Panels Sequencing kits

A seguire: **Rinfresco offerto da Diatech Labline**

Per iscriverti all'evento clicca il link o usa il QR code!

www.labline.it/evento180619



Diatech Labline srl

Via Ignazio Silone, 1/B - 60035 Jesi (AN)

Phone: 0731 213243

Email: commerciale@labline.it

